BlOlogical[™] Mail-In Spore Testing Solutions



Automatic Test Tracking with Minimal Fffort

- Steam cycle test results in 24 hours*
- Test results automatically recorded online
- Test failures are immediately communicated via e-mail or phone



Each test packet includes 1 test strip and 1 control strip Test packets are compatible with steam, dry heat, chemical and ETO. *Results in 24 hours once envelope is received at the test lab.

BIOlogical and OnTraQ

A Powerful Combination in Infection Control Compliance

- 24/7 online access to spore test results
- Automatic email and text alerts
- Practice task management and reporting

Seamless integration with HealthFirst products











Unused w Recovery



Medical Waste Recovery



Amalgam





BIOlogical^{**}10 10 Hour In-Office Spore Testing System



BIOlogical 10 Steam Biological indicators

100 indicators/Box Item# 3703

BIOlogical In-Office Monitoring Incubator

100 indicators/Box Item# 3703

A Breakthrough For In-Office **Biological Monitoring**

- Spore growth detected as few as 3-5 hours
- Final test results in 10 hours
- Access to online documentation of test results

Easy to use & Easy to read results:

Simply process the indicator in a normal sterilizer cycle. Activate the indicator, place it in the incubator, then see and record the results 10 hours later! If any bacterial spores survive the process, the indicator media will turn yellow.

Complete:

Everything you need for effective biological monitoring of steam sterilizers, including incubator, indicators, manual and access to the OnTraQ online compliance management software.

BIOlogical and OnTraQ

A Powerful Combination in Infection Control Compliance

- 24/7 online access to spore test results
- Automatic email and text alerts
- Practice task management and reporting

Seamless integration with HealthFirst products









Unused Meds Recovery



Medical Waste Recovery



Amalgam



